

Product datasheet for **SM1026F**

Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]
Clone Name:	HP6050
Applications:	ELISA, IF, IHC
Recommended Dilution:	FLISA: 1/200-1/400 Immunocytochemistry. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS Label: FITC State: Liquid IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)