

Product datasheet for **SM1026HRP**

Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]
Clone Name:	HP6050
Applications:	ELISA, FC, IF, IHC, WB
Recommended Dilution:	ELISA: 1/4,000-1/10,000. Western blotting. Dot- and Slot-immunoblotting. Immunohistochemistry on Frozen Sections: 1/10-1/20.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgG ₃ myeloma protein
Isotype:	IgG1
Formulation:	Stock solution in 50% PBS, pH 7.4, 50% Glycerol Label: HRP State: Liquid IgG fraction Preservative: None Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)