

Product datasheet for **SM1026R**

Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: HP6050]
Clone Name:	HP6050
Applications:	ELISA, IF, IHC
Recommended Dilution:	FLISA: 1 µg/ml. Immunocytochemistry. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS and a stabilizing agent Label: PE State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: R-phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Note:	To insure lot-to-lot consistency, each batch of product is tested by ELISA or FLISA for conformance to characteristics of a standard reference reagent.



[View online »](#)