

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for SM1027B

## Rat IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARA-1]

## **Product data:**

Product Type:	Secondary Antibodies
Product Name:	Rat IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARA-1]
Clone Name:	MARA-1
Applications:	ELISA
Recommended Dilution:	ELISA: 500 ng/ml.
Reactivity:	Rat
Host:	Mouse
lsotype:	lgG1
Formulation:	PBS containing 0.09% Sodium Azide and 50% Glycerol. Label: Biotin
	State: Liquid purified lgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US