

Product datasheet for **SM1027F**

Rat IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARA-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgA (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARA-1]
Clone Name:	MARA-1
Applications:	FC
Recommended Dilution:	Flow cytometry (5 µg/ml): use 50 µl of the suggested working dilution to label 10e6 cells in 100 µl.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Purified IR1060 IgA rat myeloma protein
Isotype:	IgG1
Formulation:	PBS containing 50% glycerol and 0.09% sodium azide Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)