

Product datasheet for **SM1029B**

Rat IgE Mouse Monoclonal Antibody [Clone ID: MARE-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgE Mouse Monoclonal Antibody [Clone ID: MARE-1]
Clone Name:	MARE-1
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 5 µg/ml as detecting antibody. Western Blot.
Reactivity:	Rat
Host:	Mouse
Immunogen:	IR162, IR1016, IR2 and IR410 rat IgE myeloma proteins. Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide and 50% Glycerol Label: Biotin State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on IR2 rat IgE myeloma
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)