

Product datasheet for **SM1029HRP**

Rat IgE Mouse Monoclonal Antibody [Clone ID: MARE-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgE Mouse Monoclonal Antibody [Clone ID: MARE-1]
Clone Name:	MARE-1
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 0.5 µg/ml. Western Blot.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgE myeloma protein (IR162, IR410). Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.01% Thiomersal and 50% Glycerol. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)