

## Product datasheet for SM1030B

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

Clone Name: MARG1.2

**Applications:** ELISA

Recommended Dilution: ELISA: 500ng/ml.

**Reactivity:** Rat

Host: Mouse

Immunogen: Rat IR27 IgG1 myeloma protein. Spleen cells from immunised BALB/c mice were fused with

cells of the mouse SP2/0 myeloma cell line.

Isotype: IgG1

**Formulation:** PBS, pH7.4 containing 0.09% Sodium Azide and 50% Glycerol

Label: Biotin

State: Liquid purified IgG

**Concentration:** lot specific

**Purification:** Affinity chromatography

Conjugation: Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.