

Product datasheet for SM1030F

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

Clone Name: MARG1.2

Applications: FC

Recommended Dilution: Flow Cytometry: Use 50 μl of 5 μg/ml diluted antibody to label 16 cells in 100 μl.

Reactivity: Rat

Host: Mouse

Immunogen: Rat IR27 IgG1 myeloma.

Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma

cell line.

Isotype: IgG1

Formulation: PBS containing 0.09% Sodium Azide and 50% Glycerol

Label: FITC

State: Liquid purified IgG fraction

Label: Fluorescein Isothiocyanate Isomer 1

Concentration: lot specific

Purification: Affinity Chromatography

Conjugation: FITC

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

