

Product datasheet for SM1030HRP

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG1.2]

Clone Name: MARG1.2

Applications: ELISA

 $\label{eq:Recommended Dilution: ELISA: 0.5 $\mu g/ml.$} Recommended Dilution: ELISA: 0.5 $\mu g/ml.$$

Reactivity: Rat

Host: Mouse Isotype: IgG1

Formulation: PBS, pH 7.4 containing 0.01% Thiomersal as preservative and 50% Glycerol as stabilizer.

Label: HRP

State: Liquid purified IgG fraction. Label: Horseradish Peroxidase

Concentration: lot specific

Purification: Affinity Chromatography.

Conjugation: HRP

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

