

Product datasheet for SM1031B

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]

Clone Name: MARG2a-1

Applications: ELISA

Recommended Dilution: ELISA: 500ng/ml.

Reactivity: Rat

Host: Mouse Isotype: IgG1

Formulation: PBS, pH7.4 containing 0.09% Sodium Azide and 50% Glycerol

Label: Biotin

State: Liquid purified IgG

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

