

Product datasheet for **SM1031F**

Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]
Clone Name:	MARG2a-1
Applications:	FC
Recommend Dilution:	Flow Cytometry: Use 10 µl of 5 µg/ml diluted antibody to label 10 ⁶ cells in 100 µl.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Purified Rat IgG2a Spleen cells from immunised BALB/c mice were fused with cells of the PAI-0 mouse myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide and 50% Glycerol Label: FITC State: Liquid purified IgG fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)