

## Product datasheet for **SM1031HRP**

### Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]
Clone Name:	MARG2a-1
Applications:	ELISA
Recommended Dilution:	ELISA: 0.5µg/ml.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS, pH7.4 containing 0.01% Thiomersal and 50% Glycerol Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)