

## Product datasheet for **SM1031P**

### Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2a-1]
Clone Name:	MARG2a-1
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml; This antibody may be used as a coating antibody in ELISA in combination with clone MARK-1/MARL-15 as detection reagent and purified rat IgA (SA021).
Reactivity:	Rat
Host:	Mouse
Immunogen:	Polyclonal rat IgG2a. Spleen cells from immunised BALB/c mice were fused with cells of the mouse PAI-O myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)