

Product datasheet for **SM1032F**

Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-3]
Clone Name:	MARG2b-3
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: 5 µg/ml, use 50 µl of the suggested working dilution to label 10e6 cells in 100µl.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide and 50% Glycerol Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)