

Product datasheet for **SM1033B**

Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4]

Product data:

| | |
|-----------------------|--|
| Product Type: | Secondary Antibodies |
| Product Name: | Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4] |
| Clone Name: | MARM-4 |
| Applications: | ELISA |
| Recommended Dilution: | ELISA: 125ng/ml as detecting antibody. |
| Reactivity: | Rat |
| Host: | Mouse |
| Immunogen: | IR202 rat IgM myeloma protein. Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line. |
| Isotype: | IgG1 |
| Formulation: | Phosphate buffered saline, 0.09% Sodium Azide, 50% Glycerol Label: Biotin State: Liquid purified IgG |
| Concentration: | lot specific |
| Purification: | Affinity chromatography |
| Conjugation: | Biotin |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)