

Product datasheet for SM1033HRP

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4]

Clone Name: MARM-4
Applications: ELISA

Recommended Dilution: ELISA: 0.5 µg/ml.

Reactivity: Rat

Host: Mouse

Immunogen: IR202 rat IgM myeloma protein.

Isotype: IgG1

Formulation: PBS, pH7.4 containing 0.01% Thiomersal and 50% Glycerol

Label: HRP

State: Liquid purified IgG Label: Horseradish Peroxidase

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: HRP

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.