

Product datasheet for **SM1033P**

Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARM-4]
Clone Name:	MARM-4
Applications:	ELISA
Recommended Dilution:	ELISA: 5 µg/ml as Coating antibody in combination with clone MARK-1/MARL-15 (AM05512HR-N) as Detection reagent and purified IgM (SA025).
Reactivity:	Rat
Host:	Mouse
Immunogen:	IR202 rat IgM myeloma protein. Spleen cells from immunised BALB/c mice were fused with cells of the mouse Sp2/0 myeloma cell line.
Isotype:	IgG1
Formulation:	State: Purified State: Liquid purified Ig fraction containing 0.09% Sodium Azide as preservative.
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)