

## Product datasheet for **SM1046**

### Feline IgG Mouse Monoclonal Antibody [Clone ID: K62 2G4]

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Secondary Antibodies  |
| Product Name:         | Feline IgG Mouse Monoclonal Antibody [Clone ID: K62 2G4]  |
| Clone Name:           | K62 2G4   |
| Applications:         | ELISA   |
| Recommended Dilution: | ELISA: 1/300 (detection).   |
| Reactivity:           | Feline  |
| Host:                 | Mouse   |
| Immunogen:            | Feline IgG. Spleen cells from immunised mice were fused with cells of the mouse P3-X63-Ag8.653 myeloma cell line.                                   |
| Isotype:              | IgG1  |
| Formulation:          | containing 0.2M Tris/HCl pH7.4 and 5-10% foetal calf serum and 0.01% Sodium Azide<br>State: Supernatant<br>State: Liquid Tissue Culture Supernatant |
| Conjugation:          | Unconjugated  |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.                   |



[View online »](#)