

Product datasheet for SM1349P

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Actin (F-Actin) Mouse Monoclonal Antibody [Clone ID: NH3]

Product data:

Product Type: Primary Antibodies

Clone Name: NH3

Applications: ELISA, FC, IF, IHC, WB

Recommended Dilution: ELISA: 1/10.

Western Blotting: 1/100-1/500.

Flow Cytometry: 1/10 Immunofluorescence.

Immunohistochemistry on Frozen Sections.

Reactivity: Human, Mouse, Rabbit, Rat

Host: Mouse Isotype: IgM

Clonality: Monoclonal

Immunogen: Human monocytes and U937 cell line.

Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma

cell line.

Specificity: This antibody recognizes Human Filamentous Actin (F-actin).

The antibody binds to the N-terminal region of Actin, butnot to the extreme N-terminal 40 amino acids. In tissue sections the antibody stains the cytoplasm of macrophages strongly,

and gives granular, localised nuclear staining of all cell types.

Clone NH3 is reported to recognize Actin in the filamentous form with the epitope likely to be

located between residues 120 and 226 of the molecule.

Clone NH3 is also described to show reactivity with a 43 kDa polypeptide using cell lines U937

and HL-60 by SDS/PAGE and Immunoblotting.

Formulation: PBS

State: Purified

State: Liquid purified IgM fraction prepared from Tissue Culture Supernatant

Preservative: 0.09% Sodium Azide

Concentration: lot specific

Conjugation: Unconjugated





Actin (F-Actin) Mouse Monoclonal Antibody [Clone ID: NH3] - SM1349P

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Background: Actins are highly conserved proteins that are involved in various types of cell motility and are

ubiquitously expressed in all eukaryotic cells. G-Actin (globular actin) with bound ATP can polymerise to form F-Actin (Filamentious Actin). Actin monomers spiral ariound the axis of

the filament similar to a double helix.

F-Actin may also undergo a process called treadmilling, in which filament length remains

constant and actin monomers add at one end and dissociated at the other.

Synonyms: Actin F type