

## Product datasheet for **SM1490B**

### Mouse IgG Rat Monoclonal Antibody [Clone ID: cocktail]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG Rat Monoclonal Antibody [Clone ID: cocktail]
Clone Name:	cocktail
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1 µg/ml.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG from BALB/c mice. Spleen cells from immunised LOU/c rats were fused with cells of the IR983F rat myeloma cell line.
Isotype:	cocktail
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction from Tissue Culture Supernatant Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Mouse Ig/Sepharose 4B
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)