

## Product datasheet for **SM1491P**

### Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-13]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG1 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG1-13]
Clone Name:	LO-MG1-13
Applications:	ELISA
Recommended Dilution:	ELISA: 5 µg/ml as coating antibody. This antibody may be used as a Coating antibody in ELISA in combination with SM1017HRP as the detection reagent.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified mouse IgG1 from BALB/c mice
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)