

## Product datasheet for **SM1492HRP**

### Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-9]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-9]
Clone Name:	LO-MG2a-9
Applications:	ELISA
Recommended Dilution:	ELISA (250 ng/ml).
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified Mouse IgG2a heavy chain.
Isotype:	IgG1
Formulation:	PBS containing 50% glycerol as stabilizer and 0.01% Thimerosal as preservative Label: HRP State: Liquid IgG fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)