

Product datasheet for **SM1492P**

Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-9]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG2a-9]
Clone Name:	LO-MG2a-9
Applications:	ELISA
Recommended Dilution:	ELISA: 5µg/ml; This antibody may be used as a coating antibody in ELISA in combination with LO-MG2a-7 (SM1018HRP) as the detection reagent.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Purified IgG2a from BALB/c mice. Spleen cells from immunised LOU/c rat were fused with cells of the IR983F rat myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)