

Product datasheet for **SM1494P**

Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-13]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG3 (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MG3-13]
Clone Name:	LO-MG3-13
Applications:	ELISA
Recommended Dilution:	ELISA (5 µg/ml): This antibody may be used as a Coating antibody in ELISA in combination with LO-MG3-7 (SM1022B) as the Detection reagents.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Polyclonal mouse IgG. Spleen cells from immunised Lou/c rats were fused with cells of the rat IR983F myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on mouse IgG3/Sepharose 4B
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)