

Product datasheet for **SM1495P**

Mouse IgM (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MM-3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM (heavy chain) Rat Monoclonal Antibody [Clone ID: LO-MM-3]
Clone Name:	LO-MM-3
Applications:	ELISA
Recommended Dilution:	ELISA: 5 µg/ml (Coating). This antibody may be used as a Coating antibody in Sandwich ELISA in combination with LO-MM-9 (SM1025B) and Streptavidin-HRP as the detection reagents.
Reactivity:	Mouse
Host:	Rat
Immunogen:	Polyclonal Mouse IgM.
Isotype:	IgM
Formulation:	PBS containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgM fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Mouse IgM/Sepharose 4B.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)