

Product datasheet for **SM1496**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: LO-MK-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	LO-MK-1
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Useful as second antibody in Immunoassays, such as ELISA (good binding on plastic), Immunohistochemistry, Flow Cytometry, Immunofluorescence and Western blotting.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognizes a Mouse Gamma Light Chain of Immunoglobulin (determined by Immunodot) and exhibits no reactivity to Mouse Lambda Light Chain. No detectable reactivity with Human IgG as tested by ELISA.
Formulation:	State: Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.



[View online »](#)