

Product datasheet for **SM1496F**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: LO-MK-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	LO-MK-1
Applications:	ELISA, IHC
Recommended Dilution:	Useful as second antibody in Immunoassays , such as ELISA (good binding on plastic), Immunohistochemistry on Frozen Sections and Flow Cytometry .
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognizes a Mouse Gamma Light Chain of Immunoglobulin (determined by Immunodot) and exhibits no reactivity to Mouse Lambda Light Chain. No detectable reactivity with Human IgG as tested by ELISA.
Formulation:	PBS Label: FITC State: Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)