

Product datasheet for **SM1496P**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: LO-MK-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	LO-MK-1
Applications:	ELISA, IHC
Recommended Dilution:	ELISA: 500-1000 ng/ml. Immunohistochemistry on Frozen Sections.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Purified Mouse IgG1 kappa
Specificity:	This antibody recognizes Murine Kappa Immunoglobulin Light Chain, and does not cross-react with other Murine Immunoglobulin classes or subclasses.
Formulation:	PBS Label: HRP State: Purified State: Liquid purified IgG fraction Stabilizer: 50% glycerol Preservative: 0.01% Thiomersal Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)