

Product datasheet for **SM1497B**

Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-8]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-8]
Clone Name:	MARG2b-8
Applications:	ELISA
Recommended Dilution:	ELISA: 500ng/ml.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgG2b myeloma protein (IR863). Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Isotype:	IgG1
Formulation:	PBS, pH7.4 containing 0.09% Sodium Azide and 50% Glycerol Label: Biotin State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)