

Product datasheet for **SM1497P**

Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-8]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (heavy chain) Mouse Monoclonal Antibody [Clone ID: MARG2b-8]
Clone Name:	MARG2b-8
Applications:	ELISA
Recommended Dilution:	ELISA: 5 µg/ml (as a coating antibody).
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgG2b myeloma protein (IR863). Spleen cells from immunised BALB/c mice were fused with cells of the mouse Sp2/0 myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)