

Product datasheet for **SM1498P**

Canine IgE Mouse Monoclonal Antibody [Clone ID: E6-2A1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgE Mouse Monoclonal Antibody [Clone ID: E6-2A1]
Clone Name:	E6-2A1
Applications:	ELISA, WB
Recommended Dilution:	Western Blot. ELISA: The antibody may be used as a Capture antibody in ELISA in combination with SM1499B as detecting antibody.
Reactivity:	Canine
Host:	Mouse
Immunogen:	Affinity purified IgE preparation from heavily parasitized and allergic Canine (Dog) Serum.
Isotype:	IgG2a
Formulation:	PBS, pH 7.2 containing 0.05% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)