

Product datasheet for **SM1499P**

Canine IgE Mouse Monoclonal Antibody [Clone ID: E6-71A1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgE Mouse Monoclonal Antibody [Clone ID: E6-71A1]
Clone Name:	E6-71A1
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blotting: Under reducing conditions, this antibody detects a single major band of approximately 62-65 kDa in samples of affinity purified Dog IgE.
Reactivity:	Canine
Host:	Mouse
Immunogen:	Affinity purified IgE preparation from heavily parasitized and allergic dog serum.
Isotype:	IgG1
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgG fraction Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A from tissue culture supernatant
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)