

Product datasheet for **SM1500**

Guinea pig IgG Mouse Monoclonal Antibody [Clone ID: MsGp3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Guinea pig IgG Mouse Monoclonal Antibody [Clone ID: MsGp3]
Clone Name:	MsGp3
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections.
Reactivity:	Guinea Pig
Host:	Mouse
Immunogen:	In vivo passaged EN-L2C guinea pig leukaemia B cell line. Spleen cells from an immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
Isotype:	IgM
Formulation:	containing 0.2M Tris/HCl pH7.2 and 0.09% Sodium Azide State: Supernatant State: Liquid Tissue Culture Supernatant
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)