

Product datasheet for **SM1508P**

Feline IgA Mouse Monoclonal Antibody [Clone ID: IgA9-6A]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgA Mouse Monoclonal Antibody [Clone ID: IgA9-6A]
Clone Name:	IgA9-6A
Applications:	ELISA, WB
Recommended Dilution:	Western Blot. ELISA: This antibody may be used as a capture antibody in an ELISA in combination with either SM1503B or AM01151BT-N as detecting antibody.
Reactivity:	Feline
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS, pH 7.2 containing 0.05% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)