

Product datasheet for SM1526F

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD81 / TAPA1 Hamster Monoclonal Antibody [Clone ID: Eat2]

Product data:

Product Type: Primary Antibodies

Clone Name: Eat2
Applications: FC

Recommended Dilution: Flow Cytometry: The Fc region of monoclonal antibodies may bind non-specifically to cells

expressing low affinity Fc receptors.

Reactivity: Mouse
Host: Hamster
Isotype: IgG

Clonality: Monoclonal

Immunogen: 38C13, murine B cell line. Spleen cells from immunised Armenian hamsters were fused with

cells of the mouse PX3-Ag.8.653 myeloma cell line.

Specificity: This antibody recognises and CD81. Eat 2 induces homotypic aggregation of B cells and

inhibits anti Ig and IL-4 induced proliferation. Eat 2 requires the presence of both extracellular

loops of TAPA-1 for binding.

Formulation: PBS, pH 7.4, containing 1% BSA as stabilizer and 0.09% Sodium Azide as preservative.

Label: FITC

State: Liquid purified IgG fraction.

Label: Fluorescein Isothiocyanate Isomer 1

Concentration: lot specific

Purification: Affinity Chromatography on Protein G.

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: CD81 antigen

Database Link: Entrez Gene 12520 Mouse

P35762






Background:

CD81 is a 26kD cell surface glycoprotein that is also known as TAPA-1. In rodents CD81 is expressed at much higher levels on resting B cells than on T cells, although increased expression on T cells is found following activation.

Synonyms:

S5.7; TAPA-1; TAPA1; Tetraspanin-28; Tspan-28; TSPAN28

Product images:



Staining of Mouse spleen cells with Hamster Anti Mouse CD81-FITC (SM1526F/FT).