

## **Product datasheet for SM1602P**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### CD147 (BSG) Mouse Monoclonal Antibody [Clone ID: MEM-M6/1]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: MEM-M6/1

**Applications:** FC, IF, IHC, IP, WB

Recommended Dilution: Flow Cytometry: Use 10 μl of 1/50-1/100 diluted antibody to label 16 cells or 100 μl whole

blood.

Immunoprecipitation.

**Western Blot** (under non-reducing conditions).

Immunohistochemistry on Paraffin Sections: 10 µg/ml.

This Basigin antibody was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody, followed by alkaline phosphatase-streptavidin and

chromogen.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: CD147Rg, consisting of the cDNA coding for the hinge region, CH2 and CH3 domains of

Human IgG1.

**Specificity:** This antibody recognises the CD147 cell surface antigen.

**Formulation:** PBS, pH 7.4

State: Purified

State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide

**Concentration:** lot specific

**Purification:** Affinity Chromatography on Protein G

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.





### CD147 (BSG) Mouse Monoclonal Antibody [Clone ID: MEM-M6/1] - SM1602P

**Stability:** Shelf life: one year from despatch.

Gene Name: basigin (Ok blood group)

Database Link: Entrez Gene 682 Human

P35613


**Background:** CD147 is a 50-60kD glycoprotein also known as neurothelin and basigin. CD147 is broadly

expressed on peripheral blood cells and by endothelial cells. Expression upon T cells is


increased following activation.

Synonyms: Basigin, BSG, EMMPRIN, TCSF

# **Product images:**




Immunohistochemistry: CD147 Antibody staining Formalin-Fixed, Paraffin-Embedded Human Placenta at 10 ug/ml.



Immunohistochemistry: CD147 Antibody staining Formalin-Fixed, Paraffin-Embedded Human Prostate at 10 ug/ml.





Flow Cytometry: Staining of Human Erythrocytes with Mouse anti Human CD147 Antibody