

Product datasheet for **SM1814PSX**

Human IgE Mouse Monoclonal Antibody [Clone ID: 4C3]


Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 4C3]
Clone Name:	4C3
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1 µg/ml. Sandwich ELISA: 5 µg/ml for coating plates. Can be used as a Captive antibody for Sandwich ELISA (paired with 1A2 HRP conjugated). Western Blot: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Immunogen:	Recombinant Human IgE
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Aff - Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)

Product images:



Sandwich ELISA using 4C3 coated plate and HRP conjugated 1A2