

## Product datasheet for **SM2066R**

### CD4 Mouse Monoclonal Antibody [Clone ID: 2-35]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	2-35
Applications:	FC
Recommended Dilution:	Flow Cytometry: Neat - 1/5; Use 10µl of the suggested working dilution to label 10e6 cells in 100µl.
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Chicken embryonic thymocytes. Lymph node cells from immunised BALB/c mice were fused with cells of the SP2/0 myeloma cell line.
Specificity:	This antibody recognises the homologue of human CD4, a 64kD cell surface protein expressed by thymocytes and a subset of T cells.
Formulation:	PBS, pH 7.4 containing 0.09% Sodium Azide and 1% Bovine Serum Albumin Label: PE State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	PE
Storage:	Prior to and following reconstitution store the antibody at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">Q6Q4G0</a>



[View online »](#)

**Background:** CD4 is a single chain transmembrane glycoprotein (59 kDa) which belongs to the immunoglobulin superfamily. CD4 is present on a subset of T lymphocytes ("helper/inducer" T cells) and is also expressed at a lower level on monocytes, tissue macrophages and granulocytes. The antigen is involved in binding to MHC class II molecules. The intracellular domain of the antigen is associated with p56lck protein tyrosine kinase.

**Synonyms:** T-cell surface antigen T4/Leu-3