

Product datasheet for **SM2069AP**

IGHG4 Mouse Monoclonal Antibody [Clone ID: HP6025]

Product data:

Product Type:	Secondary Antibodies
Product Name:	IGHG4 Mouse Monoclonal Antibody [Clone ID: HP6025]
Clone Name:	HP6025
Applications:	ELISA, FC, IF, IHC, IP, WB
Recommended Dilution:	ELISA: 1/1000-1/2000 (2-9). Immunohistochemistry on Frozen Sections: reported in literature (5,7,10,11). Immunohistochemistry on Paraffin Sections: reported in literature (10). Flow Cytometry: reported in literature (15, 18, 19). Immunocytochemistry: reported in literature (12). Western Blot: reported in literature (3,13-16).
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG4 myeloma protein
Isotype:	IgG1
Formulation:	50mM Tris/ 1 mM MgCl ₂ , 50 % Glycerol, pH 8.0 Label: AP State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Conjugation:	AP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)