

Product datasheet for SM2094P

Bsg Rat Monoclonal Antibody [Clone ID: OX114]

Product data:

Product Type: Primary Antibodies Clone Name: OX114 **Applications:** FC, IHC, WB Recommended Dilution: Western Blot. Flow Cytometry: Use 10 µl of 1/25-1/50 diluted antibody to label 10e6 cells in 100 µl. This antibody is routinely tested in flow cytometry on Mouse splenocytes. Immunohistochemistry on Frozen Sections. **Reactivity:** Mouse Host: Rat Isotype: lgG1 Monoclonal **Clonality:** Immunogen: Fusion protein mCD147dO-2rCD4d3 +4. Specificity: This antibody recognises CD147 (also known as basigin). OX-114 recognises both isoforms of CD147, the common form containing 2 Ig-like domains and the rare form containing 3 Ig-like domains. Note that for staining of whole blood samples increased amounts of antibody may be required due to the presence of CD147 on erythrocytes. Formulation: PBS, pH7.4 containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction. **Concentration:** lot specific **Purification:** Affinity Chromatography on Protein G. **Conjugation:** Unconjugated Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Storage: Avoid repeated freezing and thawing. Stability: Shelf life: one year from despatch. Gene Name: basigin Database Link: Entrez Gene 12215 Mouse P18572



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Bsg Rat Monoclonal Antibody [Clone ID: OX114] - SM2094P

Background:CD147 is a cell surface glycoprotein that is widely expressed on leukocytes, red cells, platelets
and endothelial cells. The reported functions of CD147 are diverse and include involvement in
T-cell activation, cell adhesion and induction of matrix metalloproteinases.

Synonyms: Basigin, BSG, EMMPRIN, TCSF

Product images:



Staining of Mouse peripheral blood lymphocytes with Rat anti-Mouse CD147 antibody.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US