

Product datasheet for **SM2108P**

Human IgG2 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: HP6200]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG2 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: HP6200]
Clone Name:	HP6200
Applications:	ELISA, IHC, R, WB
Recommended Dilution:	ELISA. Western blotting. Radioimmunoassay. Immunohistochemistry on Frozen Sections: The epitope recognized by this antibody is reported to be sensitive to formaldehyde fixation and tissue processing. We recommend the use of acetone fixation for frozen sections. Does not work on Paraffin Sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Semi-purified normal Human IgG2
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified IgG fraction from Tissue Culture Supernatant Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)