

Product datasheet for **SM2203P**

Rabbit IgM (B Cell Marker) Mouse Monoclonal Antibody [Clone ID: NRBM]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgM (B Cell Marker) Mouse Monoclonal Antibody [Clone ID: NRBM]
Clone Name:	NRBM
Applications:	FC, IHC
Recommended Dilution:	Immunohistochemistry on Frozen Sections. Flow Cytometry: Use 10 µl of diluted antibody to label 1x10 ⁶ cells in 100 µl.
Reactivity:	Rabbit
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G from tissue culture supernatant
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Should this product contain a precipitate we recommend microcentrifugation before use.



[View online »](#)