

Product datasheet for **SM2264PS**

Dpp4 Rat Monoclonal Antibody [Clone ID: H207-1082]

Product data:

Product Type:	Primary Antibodies
Clone Name:	H207-1082
Applications:	ELISA, FC
Recommended Dilution:	Flow Cytometry: Use 10 µl of 1/25-1/100 diluted antibody to label 106 cells in 100 µl. ELISA.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognises CD26, which is also known as thymocyte-activating molecule (THAM).
Formulation:	PBS, pH 7.4 containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	dipeptidylpeptidase 4
Database Link:	Entrez Gene 13482 Mouse P28843



[View online »](#)

Background:

Murine CD26 is a type II integral membrane dipeptidyl-peptidase that is expressed at the cell surface as a heterodimer, of approximately 220 kDa, and has a soluble form in the plasma. The CD26 molecule is a multifunctional glycoprotein that has ectopeptidase activity and is also reported to be involved in cell adhesion, signal transduction in T lymphocytes and HIV infection. CD26 is differentially expressed on developing thymocytes and is upregulated on activated T-cells but is not expressed on resting and thioglycollate stimulated macrophages. Expression of CD26 has also been reported on cells from a variety of non-haematopoietic tissues including liver, intestine and kidney.

Synonyms:

ADCP2, Dipeptidyl peptidase 4, DPP IV, Dipeptidyl peptidase IV, ADABP, TP103