

## Product datasheet for **SM2273P**

### Chicken IgG / Chicken IgY (heavy chain) Mouse Monoclonal Antibody [Clone ID: 16C7]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (heavy chain) Mouse Monoclonal Antibody [Clone ID: 16C7]
Clone Name:	16C7
Applications:	ELISA, FC, WB
Recommended Dilution:	Flow Cytometry. ELISA. Western Blot under non-reducing conditions.
Reactivity:	Duck
Host:	Mouse
Immunogen:	Purified Pekin duck yolk IgY. Spleen cells from immunised Balb/c mice were fused with cells of the SP2/0 mouse myeloma cell line.
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)