

Product datasheet for **SM2292P**

Yellow Fever Virus Mouse Monoclonal Antibody [Clone ID: 3576]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3576
Applications:	ELISA
Recommended Dilution:	ELISA: 1/20 - 1/200.
Reactivity:	Yellow Fever Virus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognises the Yellow Fever Virus. It is negative for other Flaviviruses tested as shown by ELISA.
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)