

## Product datasheet for **SM273PS**

### CD44 Mouse Monoclonal Antibody [Clone ID: OX-50]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OX-50
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: 1/50; Use 10µl of the suggested working dilution to label 10e6 cells in 100µl. Immunohistochemistry on frozen sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Rat T cell blasts. Spleen cells from immunised BALB/c mice were fused with cells of the mouse NSO/1 myeloma cell line.
Specificity:	This antibody recognises the rat CD44 cell surface antigen.
Formulation:	PBS, pH7.4 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">P26051</a>
Background:	CD44 is expressed by T cells, B cells, macrophages and thymocytes. CD44 is a glycoprotein of 85kD, with expression being increased following activation.
Synonyms:	LHR, MDU2, MDU3, MIC4, CDw44, Epican, ECMR-III, HUTCH-I, Heparan sulfate proteoglycan, Hermes antigen, Hyaluronate receptor, PGP-1



[View online »](#)