

Product datasheet for SM3023R

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

LOC105369261 Mouse Monoclonal Antibody [Clone ID: MEM-59]

Product data:

Product Type: Primary Antibodies

Clone Name: MEM-59
Applications: FC, IHC, WB

Recommended Dilution: Flow Cytometry analysis of human blood cells using 20 μl reagent / 100 μl of whole blood or

10e6 cells in a suspension.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human T lymphocytes

Specificity: This antibody recognizes neuraminidase-sensitive epitope on CD43 (Leukosialin), a 95-135

kDa type I transmembrane glycoprotein (mucin-type) which is involved in lymphocyte

activation. CD43 is expressed by platelets and at high levels on the surface of all leukocytes; it

is negative on resting B lymphocytes and erythrocytes.

Formulation: Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade

protease free Bovine Serum Albumin (BSA) as a stabilizing agent

Label: PE

State: Liquid Ig fraction

Label: Conjugated with R-Phycoerythrin

Purification: Size-exclusion chromatography

Conjugation: PE

Storage: Store the antibody at 2 - 8 °C. DO NOT FREEZE! This product is photosensitive and should be

protected from light.

Stability: Shelf life: one year from despatch.

Database Link: Entrez Gene 105369261 Human

P16150






Background:

CD43 (leukosialin, sialophorin) is a transmembrane mucin-like protein with high negative charge, expressed on the surface of most hematopoietic cells. CD43 contributes to a repulsive barrier that interferes with cellular adhesion, however, in certain cases also promotes leukocyte aggregation. By interaction with actin-binding proteins ezrin and moesin CD43 plays a regulatory role in remodeling T-cell morphology and regulates cell-cell interactions during lymphocyte traffic. CD43 signaling both enhances LFA-1 adhesiveness and counteracts LFA-1 induction via other receptors. Expression of CD43 causes induction of functionally active tumour suppressor p53 protein, but in case of p53 and ARF defficiency CD43 promotes tumour proliferation and viability. It appears to be an important modulator of leukocyte functions.


Synonyms:

Leukocyte sialoglycoprotein, Sialophorin, Galactoglycoprotein, SPN

Product images:



Western Blotting analysis (non-reducing conditions) of isolated peripheral blood lymphocytes of various species using anti-human CD43 (MEM-59). Lane 1: lysate of human PBL Lane 2: lysate of canine PBL Lane 3: lysate of porcine PBL



Immunohistochemistry staining of human spleen (paraffin sections) using anti-CD43 (MEM-59).