

Product datasheet for **SM3032AC**

LFA3 (CD58) Mouse Monoclonal Antibody [Clone ID: MEM-63]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MEM-63
Applications:	FC
Recommended Dilution:	Flow Cytometry analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Reactivity:	Human, Porcine
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	NALM-6 human pre-B cell line
Specificity:	This antibody reacts with CD58 (LFA-3), a 40-70 kDa glycoprotein distributed over many tissues, leukocytes, erythrocytes, endothelial cells, epithelial cells and fibroblasts.
Formulation:	Tris buffered saline (TBS) solution containing 15 mM sodium azide Label: APC State: Liquid purified Ig fraction Label: Conjugated with cross-linked Allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography and adjusted for direct use.
Conjugation:	APC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This products is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD58 molecule
Database Link:	Entrez Gene 965 Human P19256




[View online »](#)

Background:

CD58 (LFA-3) is an immunoglobulin family adhesion molecule expressed by both hematopoietic and non-hematopoietic cells (often on antigen presenting cells) and serving as ligand of CD2. This interaction is important for T cell-mediated immunity. CD58 is expressed in transmembrane form and in GPI-anchored form; the later is constitutively associated with protein kinases whereas the transmembrane form activates kinase activity upon triggering. CD58 is a powerful tool for detection of minimal residual disease in acute lymphocytic leukemia, and for evaluation of liver damage related with hepatitis B.

Synonyms:

Ag3, LFA-3

Product images:


Surface staining of human peripheral blood cells with anti-CD58 (MEM-63) APC.