

## Product datasheet for **SM3041F**

### CD99 (CD99R Isoform) Mouse Monoclonal Antibody [Clone ID: MEM-131]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	MEM-131
Applications:	FC
Recommended Dilution:	Flow Cytometry analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10e6 cells in a suspension.
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	HPB-ALL human peripheral blood leukemia T-cell line
Specificity:	This antibody reacts with CD99R, an epitope restricted to a subset of CD99 molecule expressed on myeloid cells, NK cells and T lymphocytes.
Formulation:	Phosphate buffered saline (PBS) containing 15 mM sodium azide and 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA) as a stabilizing agent. Label: FITC State: Liquid purified Ig fraction Label: Conjugated with Fluorescein isothiocyanate
Conjugation:	FITC
Storage:	Store the antibody at 2 - 8 °C. DO NOT FREEZE! This product is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD99 molecule
Database Link:	<a href="#">Entrez Gene 4267 Human P14209</a>



[View online »](#)

- Background:** CD99 (E2, MIC2) is a transmembrane glycoprotein that is involved in regulation of T cell adhesive properties and programmed cell death distinct from typical apoptosis course. CD99 roles are specific to certain stages of T cell differentiation such as corticothymocytes. CD99R isoform expression is restricted in the haematopoietic system to T, NK and myeloid cells.
- Synonyms:** E2 antigen, 12E7, MIC2X, MIC2Y